

ProCare Dent 30 C - 5 l
Neutralizing agent, acidic, 5 l
with citric acid for optimum neutralization.



- Particularly suitable for transmission instruments
- Set up for reprocessing phases in Miele washer-disinfectors
- Very gentle on material
- Liquid concentrate

EAN: 4029577143466 / Material number: 12187250 / Old Material Number: 66915302EU1

Product category	
Thermal disinfectors	•
Applications	
Removal of limescale and other acid-soluble impurities	Good
Dispensing recommendation	
Dispensing recommendation per rinse	1 ml/l (0.1%)
pH value of solution	3.0
Main applications	
Dental practices	•
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Citric acid
Does not contain surfactants	•
Very good material compatibility	•
Suitable for standard instruments and handpieces	•
Dimensions and weight	
External dimensions, net height in mm	8 9/16 (217)
External dimensions, net width in mm	5 5/8 (142)
External dimensions, net depth in mm	9 13/16 (248)
External dimensions, gross height in mm	9 7/8 (250)
External dimensions, gross width in mm	7 1/2 (190)
External dimensions, gross depth in mm	5 15/16 (150)
Gross weight in kg	0.21 (6)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
CE	•
GHS / CLP compliant	•