

ProCare Dent 30 P - 5 l

Neutralizing agent, acidic, 5 l

For optimum neutralization with phosphoric acid based neutralising agent.



- Suitable for standard instruments (no transmission instruments)
- Set up for reprocessing phases in Miele washer-disinfectors
- Very gentle on material
- Liquid concentrate

EAN: 4029577143497 / Material number: 12187380 / Old Material Number: 66914302EU1

Product category	
Thermal disinfectors	•
Applications	
Removal of tarnishing, limescale, urine scale, and other acid-soluble impurities	Good
Dispensing recommendation	
Dispensing recommendation per rinse	1 ml/l (0.1%)
pH value of solution	2.3
Main applications	
Dental practices	•
Product characteristics	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Phosphoric acid
Not suitable for anodized aluminum	•
Does not contain surfactants	•
Excellent rinsing properties	•
Suitable for standard instruments	•
Not suitable for chrome-plated parts	•
Not suitable for nickel-plated parts	•
Not suitable for light and non-ferrous metals	•
Dimensions and weight	
External dimensions, net height in mm	8 9/16 (217)
External dimensions, net width in mm	5 5/8 (142)
External dimensions, net depth in mm	9 13/16 (248)
External dimensions, gross height in mm	5 1/8 (130)
External dimensions, gross width in mm	8 1/2 (215)
External dimensions, gross depth in mm	9 11/16 (245)
Gross weight in kg	0.28 (8)
Fill volume in ml	169.07 (5000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
CE	•
GHS / CLP compliant	•