

ProCare Dent 40 - 1 l [Typ 1]
Rinsing agent, 1 l
with pouring aid, for final rinsing of the load items.



- For improved and reduced drying
- Set up for reprocessing phases in Miele washer-disinfectors
- Easier dispensing thanks to precise pouring aid
- Liquid concentrate

EAN: 4029577155117 / Material number: 12187410 / Old Material Number: 66917102EU1

Product category	
Thermal disinfectors	•
Dispensing recommendation	
Dispensing recommendation per rinse	0.3 ml/l (0.03%)
pH value of solution	3.8
Small opening as pouring aid	•
Main applications	
Dental practices	•
Product characteristics	
Medium type	Rinsing agent
Aggregate form	Liquid concentrate
Ingredients	Methylchloroisothiazolinone ,Methylisothiazolinone,Phosphonate,Non-ionic tensides
Suitable for dental handpieces	•
Improves and shortens the drying process	•
Dimensions and weight	
External dimensions, net height in mm	9 1/2 (241)
External dimensions, net width in mm	4 3/16 (106)
External dimensions, net depth in mm	2 5/16 (58)
External dimensions, gross height in mm	9 11/16 (245)
External dimensions, gross width in mm	2 3/8 (60)
External dimensions, gross depth in mm	4 3/16 (105)
Gross weight in kg	0.04 (1)
Fill volume in ml	33.81 (1000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
CE	•
GHS / CLP compliant	•