

APLW 875

Installation kit

for retrofitting a DI water booster pump to a PLW 7111.



- Installed by a service technician

EAN: 4002516792444 / Material number: 12385630 /
Old Material Number: 69242005D

Product category

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PLW 7111

•

PLW 7111 [Safety]

•

Product category

Other accessories

•

Product characteristics

Material

Stainless steel, Plastic

Colour

Stainless Steel, Grey

Application area

Reprocessing of laboratory glassware

•

Dimensions and weight

External dimensions, gross height in in. (mm)

4 7/8 (123)

External dimensions, gross width in in. (mm)

12 7/16 (315)

External dimensions, gross depth in in. (mm)

19 15/16 (505)

Net weight in lb (kg)

5 (2.13)

Gross weight in lb (kg)

6 (2.45)