

APLW 891

Installation kit

Conductivity module for installing in PLW 7111.



- For determining the conductivity value in the wash water
- Installed by customer service

EAN: 4002516779414 / Material number: 12506460 / Old Material Number: 69242013D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 7111	•
PLW 7111 [Safety]	•
Product category	
Conversion kit	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of breathing masks	•
Reprocessing of compressed air unit carrier frames	•
Reprocessing of compressed air bottles	•
Dimensions and weight	
External dimensions, gross height in in. (mm)	3 3/8 (85)
External dimensions, gross width in in. (mm)	4 3/4 (120)
External dimensions, gross depth in in. (mm)	5 15/16 (150)
Net weight in lb (kg)	2 (0.49)
Gross weight in lb (kg)	2 (0.56)