

A 315/1 Module

For increased capacity and reprocessing of hollow instruments.



- Increases the injector capacity from 22 to 44
- Can be used in all A 151 lower baskets available on the market
- Compatible with all standard dental & ENT adapters & accessories
- Injector nozzles not included
- A 800 tubular filter for use in either injector manifold (not included)
- Clearances in lower basket A 151: H 210 +/- 30, W 383, D 479 mm

EAN: 4002516881186 / Material number: 12531630 / Old Material Number: 69631501D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 C	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for ENT instruments	•
Insert / module for dental instruments	•
Other inserts / other modules	•
Application area	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Option to use a central filter	•
Capacity	
ENT instruments [number]	22
Handpieces [number]	22
Siphons [number]	22
Lumened instruments [number]	22
Dimensions and weight	
External dimensions, net height in in. (mm)	4 11/16 (118)
External dimensions, net width in in. (mm)	19 3/8 (492)
External dimensions, net depth in in. (mm)	19 13/16 (503)
External dimensions, gross height in in. (mm)	6 1/2 (165)
External dimensions, gross width in in. (mm)	19 15/16 (505)
External dimensions, gross depth in in. (mm)	21 11/16 (550)
Net weight in lb (kg)	7 (2.90)
Gross weight in lb (kg)	9 (3.66)
Accessory pack	
A 838 key set [number]	1