

APLW 895

Installation kit

drain water cooling for retrofitting in the PLW 7111.



- For cooling drain water in temperature-sensitive drain
- Installed by customer service

EAN: 4002516888062 / Material number: 12616020 / Old Material Number: 69242016D

Product category	
SlimLine laboratory glassware washers	•
Machine affiliation	
PLW 7111	•
Product category	
Conversion kit	•
Product characteristics	
Material	Stainless steel, Plastic
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, gross height in in. (mm)	11 7/16 (290)
External dimensions, gross width in in. (mm)	4 1/16 (102)
External dimensions, gross depth in in. (mm)	12 13/16 (324)
Net weight in lb (kg)	3 (0.92)
Gross weight in lb (kg)	3 (1.09)