


A 101

Upper basket

Height adjustable, for the optimum loading of inserts.



- Vertical clearance 160 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm


EAN: 4002516828143 / Material number: 12645590 / Old Material Number: 69610101D

Product category	
Fresh-water dishwashers	•
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Machines for cleaning industrial parts	•
Washer-disinfectors, fire brigades, and emergency services	•
Machine affiliation	
PG 8063	•
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Upper basket	•
Carrier / basket carrier	•
Upper basket	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of different utensils	•
Reprocessing of air regulators	•
Dimensions and weight	
External dimensions, net height in in. (mm)	8 1/8 (206)
External dimensions, net width in in. (mm)	20 13/16 (528)
External dimensions, net depth in in. (mm)	20 3/4 (527)
External dimensions, gross height in in. (mm)	8 11/16 (220)
External dimensions, gross width in in. (mm)	21 5/16 (540)

A 101

Upper basket

Height adjustable, for the optimum loading of inserts.



EAN: 4002516828143 / Material number: 12645590 / Old Material Number: 69610101D

External dimensions, gross depth in in. (mm)	21 11/16 (550)
--	----------------

Net weight in lb (kg)	7 (3.03)
-----------------------	----------

Gross weight in lb (kg)	9 (3.68)
-------------------------	----------