


A 102

Upper basket

Height adjustable, for the optimum loading of inserts.



- Vertical clearance 205 +/- 30 mm, width 475 mm, depth 443 mm
- With integrated spray arm


EAN: 4002516855262 / Material number: 12645600 / Old Material Number: 69610201D

Product category	
Glassware Washers, laboratory	•
Thermal disinfectors	•
Cleaning machine, small, dental	•
Washer-disinfectors, medical	•
Machines for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Carrier / basket carrier	•
Upper basket	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of OR shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Capacity	
DIN mesh trays [number]	2
SPRI mesh trays [number]	1
Dimensions and weight	
External dimensions, net height in in. (mm)	8 1/8 (206)
External dimensions, net width in in. (mm)	20 13/16 (528)

A 102

Upper basket

Height adjustable, for the optimum loading of inserts.



EAN: 4002516855262 / Material number: 12645600 / Old Material Number: 69610201D

External dimensions, net depth in in. (mm)	20 3/4 (527)
External dimensions, gross height in in. (mm)	9 1/16 (230)
External dimensions, gross width in in. (mm)	21 5/16 (540)
External dimensions, gross depth in in. (mm)	21 1/2 (545)
Net weight in lb (kg)	7 (3.05)
Gross weight in lb (kg)	8 (3.79)