


E 167
Cutlery basket
For holding cutlery.




- Suitable for plate basket U 506
- 7 Compartments

EAN: 4002516854869 / Material number: 12663060 / Old Material Number: 67616702D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•

E 167
Cutlery basket
For holding cutlery.



EAN: 4002516854869 / Material number: 12663060 / Old Material Number: 67616702D

Product category	
Insert for cutlery	•
Product characteristics	
Material	Plastic
Colour	White
Number of compartments	7
Dimensions and weight	
External dimensions, gross height in in. (mm)	3 3/4 (95)
External dimensions, gross width in in. (mm)	6 3/4 (170)
External dimensions, gross depth in in. (mm)	21 11/16 (550)
Net weight in lb (kg)	1 (0.19)
Gross weight in lb (kg)	1 (0.34)