


A 103

Upper basket

For the optimum loading of inserts.




- Vertical clearance: 3 3/4 (95) in (mm), width 505 mm, depth 472 mm
- Particularly suitable for the reprocessing of DIN mesh trays
- Also suitable for further applications with low height items
- Can be used in combination with the A 202 mobile unit
- With integrated spray arm

EAN: 4002516855040 / Material number: 12664750 / Old Material Number: 69610302D

Product category	
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Machines for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Carrier / basket carrier	•
Upper basket	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of OR shoes	•
Reprocessing of surgical instruments	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Maximum vertical clearance in mm	3 3/4 (95)
Capacity	
DIN mesh trays [number]	2
SPRI mesh trays [number]	1
Dimensions and weight	
External dimensions, net height in in. (mm)	5 1/4 (133)
External dimensions, net width in in. (mm)	20 13/16 (528)
External dimensions, net depth in in. (mm)	20 13/16 (528)
External dimensions, gross height in in. (mm)	6 3/4 (170)
External dimensions, gross width in in. (mm)	21 1/2 (545)
External dimensions, gross depth in in. (mm)	22 1/16 (560)
Net weight in lb (kg)	9 (4.18)

A 103
Upper basket
For the optimum loading of inserts.



EAN: 4002516855040 / Material number: 12664750 / Old Material Number: 69610302D

Gross weight in lb (kg)

10 (4.64)