


A 104

Upper basket

For the optimum loading of instruments for MIS, arthroscopy and urology.



- Right side for holding 11 hollow instruments
- Max. length of instruments connected: 430 mm
- Additional Luer lock connectors for hollow instruments
- Left-hand side free for inserts (max. clearances H 310, W 240, D 455 mm)

EAN: 4002516854807 / Material number: 12664810 / Old Material Number: 69610401D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Basket for MIS instruments	•
Injector basket	•
Upper basket	•
Application area	
Reprocessing of MIS instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Capacity	
DIN mesh trays [number]	1
Lumened instruments [number]	11
Dimensions and weight	
External dimensions, net height in in. (mm)	13 1/8 (332)
External dimensions, net width in in. (mm)	20 13/16 (528)
External dimensions, net depth in in. (mm)	20 3/4 (526)
External dimensions, gross height in in. (mm)	16 9/16 (420)
External dimensions, gross width in in. (mm)	21 (533)
External dimensions, gross depth in in. (mm)	23 7/16 (594)
Net weight in lb (kg)	11 (4.86)
Gross weight in lb (kg)	15 (6.90)
Accessory pack	
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	2
E 448 silicone hose with female Luer Lock adapter [number]	2
E 452 injector nozzle 2.5 x 60 mm [number]	1
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	1
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	1