


A 105/1

Upper basket

Height adjustable, without injector nozzles, for various instruments.



- Right-hand side: space for 22 injectors
- Left-hand side free for inserts
- #Active internal drying of handpieces and contra-angle handpieces
- Clearances H 225 +/- 30, W 360, D 445 mm
- A 800 tubular filter can be used in both injector manifolds

EAN: 4002516854814 / Material number: 12664820 / Old Material Number: 69610502D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Injector basket	•
Upper basket	•
Application area	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Option to use a central filter	•
Capacity	
Handpieces [number]	22
Lumened instruments [number]	22
Dimensions and weight	
External dimensions, net height in in. (mm)	8 (202)
External dimensions, net width in in. (mm)	20 13/16 (528)
External dimensions, net depth in in. (mm)	20 13/16 (528)
External dimensions, gross height in in. (mm)	8 11/16 (220)
External dimensions, gross width in in. (mm)	21 5/16 (540)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	11 (4.88)
Gross weight in lb (kg)	12 (5.57)
Accessory pack	
A 838 key set [number]	1