

A 846

Injector nozzle

With support, for cleaning laboratory flasks and round flasks.



- For laboratory bottles up to 10 L and round flasks up to 6 L
- Nozzle length: 127 mm, nozzle length with holder: 158 mm
- Nozzle Ø: 10 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516855088 / Material number: 12665100 / Old Material Number: 69684601D

Product category	
Laboratory glasswashers	•
SlimLine laboratory glassware washers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
Product category	
Injector nozzle	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Rubber, Plastic
Colour	Stainless Steel, Black, Grey
Dimensions and weight	
External dimensions, net height in in. (mm)	6 1/4 (158)
External dimensions, net width in in. (mm)	4 3/8 (110)
External dimensions, net depth in in. (mm)	4 3/8 (110)
External dimensions, gross height in in. (mm)	6 5/16 (160)
External dimensions, gross width in in. (mm)	5 3/4 (145)
External dimensions, gross depth in in. (mm)	10 7/16 (265)
Net weight in lb (kg)	1 (0.23)
Gross weight in lb (kg)	1 (0.39)