

A 847

Injector nozzle

With support, for cleaning Erlenmeyer flasks.



- For Erlenmeyer flasks, up to 5 L
- Nozzle length: 127 mm
- Nozzle Ø: 10 mm, Ø of support: 142 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516855101 / Material number: 12665110 / Old Material Number: 69684701D

Product category	
Laboratory glasswashers	•
SlimLine laboratory glassware washers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
Product category	
Injector nozzle	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
Dimensions and weight	
External dimensions, net height in in. (mm)	6 3/4 (170)
External dimensions, net width in in. (mm)	9 13/16 (249)
External dimensions, net depth in in. (mm)	1 7/16 (35)
External dimensions, gross height in in. (mm)	2 13/16 (70)
External dimensions, gross width in in. (mm)	7 1/8 (180)
External dimensions, gross depth in in. (mm)	10 1/4 (260)
Net weight in lb (kg)	0 (0.13)
Gross weight in lb (kg)	1 (0.23)