

A 861  
Holding frames  
for A 301/4 and A 301/5 injector modules.



- Arrangement: 3 x 6 fields
- For use on A 301/4 and A 301/5 injector modules
- For the optimum positioning of laboratory glassware
- Part of SimpleLoad system

EAN: 4002516854876 / Material number: 12665440 / Old Material Number: 69686101D

<b>Product category</b>	
Glassware Washers, laboratory	•
SlimLine laboratory glassware washers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
<b>Machine affiliation</b>	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
<b>Product category</b>	
Other inserts / other modules	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	18
Compartment width in mm	2 15/16 (74)
Compartment length in mm	2 15/16 (74)
Part of SimpleLoad system	•
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	3 3/16 (80)
External dimensions, net width in in. (mm)	7 7/8 (200)
External dimensions, net depth in in. (mm)	13 13/16 (350)
External dimensions, gross height in in. (mm)	5 1/8 (130)
External dimensions, gross width in in. (mm)	10 1/16 (255)
External dimensions, gross depth in in. (mm)	19 1/8 (485)
Net weight in lb (kg)	1 (0.59)
Gross weight in lb (kg)	2 (0.94)