


A 862
Holding frames
for A 302/2 and A 302/3 injector modules.



- Arrangement: 4 x 8 fields
- For use on A 302/2 and A 302/3 injector modules
- For the optimum positioning of laboratory glassware
- Plastic-coated

EAN: 4002516854203 / Material number: 12665480 /
Old Material Number: 69686201D

Product category	
Glassware Washers, laboratory	•
SlimLine laboratory glassware washers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
Product category	
Other inserts / other modules	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Colour	Stainless Steel, Black
Number of compartments	32
Compartment width in mm	2 1/4 (56)
Compartment length in mm	2 1/16 (52)
Part of SimpleLoad system	•
Dimensions and weight	
External dimensions, net height in in. (mm)	1 3/16 (30)
External dimensions, net width in in. (mm)	9 7/8 (250)
External dimensions, net depth in in. (mm)	13 13/16 (350)
External dimensions, gross height in in. (mm)	3 15/16 (100)
External dimensions, gross width in in. (mm)	10 1/16 (255)
External dimensions, gross depth in in. (mm)	18 5/16 (465)
Net weight in lb (kg)	1 (0.55)
Gross weight in lb (kg)	2 (0.81)