


A 200
Mobile unit
for the optimum loading of 38 pipettes.



- for 10 100 ml pipettes in row 1, distance 20 mm
- for 14 25 ml pipettes in a 2nd row, distance 26 mm
- for 14 10 ml pipettes in a 3. row, distance 26 mm

EAN: 4002516854012 / Material number: 12667060 / Old Material Number: 69620001D

Product category	
Glassware Washers, laboratory	•
Machine affiliation	
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
Product category	
Mobile unit for pipettes	•
Application area	
Reprocessing of pipettes / phials	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Capacity	
Pipettes [number]	38
Dimensions and weight	
External dimensions, net height in in. (mm)	15 11/16 (397)
External dimensions, net width in in. (mm)	20 7/8 (529)
External dimensions, net depth in in. (mm)	21 1/2 (546)
External dimensions, gross height in in. (mm)	20 1/2 (520)
External dimensions, gross width in in. (mm)	22 1/16 (560)
External dimensions, gross depth in in. (mm)	21 11/16 (550)
Net weight in lb (kg)	11 (4.95)
Gross weight in lb (kg)	14 (6.38)