

A 202
Mobile unit
for the optimum loading of 4 DIN mesh trays.



- Suitable for 4 DIN mesh trays on 2 levels
- Insert frames: Removable for the reprocessing of tall wash items
- 6 additional connection options via Luer lock connection
- With integrated spray arm
- Loading volumes for the different levels differ from each other

EAN: 4002516854043 / Material number: 12667130 / Old Material Number: 69620201D

Product category	
Laboratory glasswashers	•
Washer-disinfectors, medical	•
Machines for cleaning industrial parts	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Capacity	
DIN mesh trays [number]	4
Dimensions and weight	
External dimensions, net height in in. (mm)	8 13/16 (223)
External dimensions, net width in in. (mm)	20 7/8 (529)
External dimensions, net depth in in. (mm)	21 3/8 (542)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	22 1/2 (570)
External dimensions, gross depth in in. (mm)	21 11/16 (550)
Net weight in lb (kg)	13 (6.05)
Gross weight in lb (kg)	15 (6.76)