


E 105/1

Insert

For the optimum loading of test tubes up to 12 x 165 mm.



- Suitable for approx. 200 Test tubes (up to 12 x 165 mm)
- Capacity including lid A 13
- The mesh width is 9 x 9 mm

EAN: 4002516853978 / Material number: 12667380 / Old Material Number: 69510503D


Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory glassware washers	•
Machines for cleaning industrial parts	•
Large-chamber laboratory glasswashers for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert / module for test tubes	•
Application area	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Capacity	
Test tubes [number]	200
Dimensions and weight	
External dimensions, net height in in. (mm)	8 3/8 (212)
External dimensions, net width in in. (mm)	7 7/8 (200)
External dimensions, net depth in in. (mm)	8 5/16 (210)
External dimensions, gross height in in. (mm)	9 5/16 (235)
External dimensions, gross width in in. (mm)	9 1/2 (240)
External dimensions, gross depth in in. (mm)	9 7/8 (250)

E 105/1

Insert

For the optimum loading of test tubes up to 12 x 165 mm.

EAN: 4002516853978 / Material number: 12667380 / Old Material Number: 69510503D



Net weight in lb (kg)	2 (0.88)
-----------------------	----------

Gross weight in lb (kg)	3 (1.16)
-------------------------	----------

<b>Accessory pack</b>	
-----------------------	--

A 13 lids [number]	1
--------------------	---