


E 139/1

Insert

For the optimum loading of test tubes up to 12 x 200 mm.



- Suitable for test tubes up to 12 x 200 mm
- Divided into 6 surfaces
- Mesh width 9 x 9 mm
- For use in lower basket
- Including lid A 13

EAN: 4002516854036 / Material number: 12667390 / Old Material Number: 69513903D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory glassware washers	•
Machines for cleaning industrial parts	•
Large-chamber laboratory glasswashers for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
Product category	
Insert / module for test tubes	•
Application area	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Capacity	
Test tubes [number]	200
Dimensions and weight	
External dimensions, net height in in. (mm)	9 5/8 (243)
External dimensions, net width in in. (mm)	7 7/8 (200)
External dimensions, net depth in in. (mm)	9 1/16 (230)
External dimensions, gross height in in. (mm)	9 11/16 (245)
External dimensions, gross width in in. (mm)	9 11/16 (245)
External dimensions, gross depth in in. (mm)	10 7/16 (265)

E 139/1

Insert

For the optimum loading of test tubes up to 12 x 200 mm.

EAN: 4002516854036 / Material number: 12667390 / Old Material Number: 69513903D



Net weight in lb (kg)	2 (1.01)
-----------------------	----------

Gross weight in lb (kg)	3 (1.21)
-------------------------	----------

<b>Accessory pack</b>	
-----------------------	--

A 13 lids [number]	1
--------------------	---