

A 203

Mobile unit

for the optimum loading of MIS instruments.



- Suitable for holding and others for E 451, E 457, E 460, E 473, E 444
- Loading on 2 levels

EAN: 4002516853985 / Material number: 12667770 / Old Material Number: 69620301D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mobile unit for MIS instruments	•
Application area	
Reprocessing of MIS instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	20 1/2 (520)
External dimensions, net width in in. (mm)	20 7/8 (529)
External dimensions, net depth in in. (mm)	21 3/8 (542)
External dimensions, gross height in in. (mm)	21 5/16 (540)
External dimensions, gross width in in. (mm)	22 1/16 (560)
External dimensions, gross depth in in. (mm)	22 1/16 (560)
Net weight in lb (kg)	20 (9.16)
Gross weight in lb (kg)	24 (11.00)
Accessory pack	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	3
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	15
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	5
E 444 drum for fibre optic cables / suction hoses [number]	1
E 445 – 12 irrigation sleeve covers E 442, opening with Ø 6 mm [number]	1
E 446 – 12 irrigation sleeve covers E 443, opening with Ø 10 mm [number]	1
E 447 female adapter, for male Luer Lock [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	6
E 449 male adapter, for female Luer Lock [number]	5
E 451 insert 1/6 mesh tray [number]	1
E 452 injector nozzle 2.5 x 60 mm [number]	3
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	8
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	6

A 203
Mobile unit
for the optimum loading of MIS instruments.

EAN: 4002516853985 / Material number: 12667770 / Old Material Number: 69620301D



E 456 spring stay for MIS instruments [number]

4

E 464 holder for injector nozzle E 454 [number]

3