

A 204  
Mobile unit  
for the optimum loading of Micro instruments.



- Level 1: suitable for inserts
- Level 2 right: 22 connections for hollow instruments
- Level 2 left: 11 Luer lock adapters, male, horizontal Hose
- With integrated spray arm
- Tubular filter A 800 can be used in both injector rails

EAN: 4002516854821 / Material number: 12667790 / Old Material Number: 69620401D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Mobile unit for ophthalmology instruments	•
<b>Application area</b>	
Reprocessing of ophthalmology instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Option to use a central filter	•
<b>Capacity</b>	
Surgery sets for micro instruments [number]	2-4
<b>Accessory pack</b>	
E 476 holder [number]	1
E 477 stopper [number]	1