

A 207

Mobile unit

For the optimum loading of Micro instruments.



- Efficient reprocessing: fast loading and high throughput
- Top level: up to 8 trays for hollow instrument loading
- Lower level: 8 horizontal trays for further ophthal. instruments

EAN: 4002516854210 / Material number: 12667800 / Old Material Number: 69620702D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mobile unit for ophthalmology instruments	•
Application area	
Reprocessing of ophthalmology instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Option to use a central filter	•
Capacity	
Trays [quantity]	16
Surgery sets for micro instruments [number]	8
Dimensions and weight	
External dimensions, net height in in. (mm)	17 1/4 (438)
External dimensions, net width in in. (mm)	20 7/8 (530)
External dimensions, net depth in in. (mm)	21 3/8 (542)
External dimensions, gross height in in. (mm)	27 1/4 (692)
External dimensions, gross width in in. (mm)	27 (685)
External dimensions, gross depth in in. (mm)	27 (685)
Net weight in lb (kg)	21 (9.58)
Gross weight in lb (kg)	29 (13.32)
Accessory pack	
A 311 tray holder [quantity]	2
A 800 tubular filter [quantity]	1
A 804 cleaning brush [quantity]	1