

A 208

Mobile unit

For the optimum loading of 8 mesh trays.



- High capacity: simultaneous reprocessing of 8 DIN mesh trays
- 4 additional connections for hollow instruments
- Removable level increases vertical clearance if required
- Information about clearances on the web: Level 1 (bottom) to level 4

EAN: 4002516854494 / Material number: 12668770 / Old Material Number: 69620801D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Mobile unit for mesh trays	•
<b>Application area</b>	
Reprocessing of mesh trays	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Capacity</b>	
DIN mesh trays [number]	8
Lumened instruments [number]	4
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	15 11/16 (397)
External dimensions, net width in in. (mm)	20 7/8 (529)
External dimensions, net depth in in. (mm)	21 1/2 (546)
External dimensions, gross height in in. (mm)	22 1/2 (570)
External dimensions, gross width in in. (mm)	18 15/16 (480)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	22 (9.83)
Gross weight in lb (kg)	24 (11.03)