

A 300/2

Injector module

For optimum loading of laboratory glassware – can be equipped with 8 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 8 injector nozzles
- Part of SimpleLoad system
- Nozzle spacing 118 mm across the width, 118 mm across the depth
- Particularly suitable for A 840 and A 841 (nozzle diameter 1/4" (6 mm))

EAN: 4002516853954 / Material number: 12668780 /  
Old Material Number: 69630004D

<b>Product category</b>	
Glassware Washers, laboratory	•
SlimLine laboratory glassware washers	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
<b>Machine affiliation</b>	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
<b>Product category</b>	
Insert / module for narrow necked glassware	•
Insert for various bottles	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel, Black
Part of SimpleLoad system	•
<b>Capacity</b>	
Narrow necked glassware [number]	8
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	2 7/8 (72)
External dimensions, net width in in. (mm)	5 1/4 (133)
External dimensions, net depth in in. (mm)	18 7/8 (479)
External dimensions, gross height in in. (mm)	7 7/8 (200)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	3 (1.25)
Gross weight in lb (kg)	4 (1.65)