


A 303
Module
for the optimum loading of up to 98 pipettes.



- Suitable for 98 measuring and full pipettes, for example
- Height of holding frame 150 mm
- 1 unit can be inserted into each of the baskets A 150 and APLW 106
- Vertical clearance: max. 450 mm (in A 150) and 605 mm (in APLW 106)

EAN: 4002516854746 / Material number: 12669220 / Old Material Number: 69630301D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Insert / module for pipettes	•
Application area	
Reprocessing of pipettes / phials	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Capacity	
Pipettes [number]	98
Dimensions and weight	
External dimensions, net height in in. (mm)	7 5/16 (185)
External dimensions, net width in in. (mm)	8 3/4 (222)
External dimensions, net depth in in. (mm)	18 9/16 (471)
External dimensions, gross height in in. (mm)	8 11/16 (220)
External dimensions, gross width in in. (mm)	10 1/4 (260)
External dimensions, gross depth in in. (mm)	19 5/16 (490)
Net weight in lb (kg)	6 (2.50)
Gross weight in lb (kg)	7 (3.08)