

## A 304 Module

For the optimum loading of up to 98 vials.



- Suitable for 98 centrifuge tubes, phials, test tubes, for example
- The loading heights differ depending on the positioning
- 1 unit can be inserted into each of the baskets A 100, A 150, and APLW 106

EAN: 4002516854982 / Material number: 12669230 /  
Old Material Number: 69630401D

<b>Product category</b>	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machines for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
<b>Machine affiliation</b>	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
PLW 7111 [IPC]	•
<b>Product category</b>	
Insert / module for test tubes	•
<b>Application area</b>	
Reprocessing of pipettes / phials	•
Reprocessing of test tubes	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Capacity</b>	
Vials [number]	98
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	5 1/8 (130)
External dimensions, net width in in. (mm)	8 3/4 (222)
External dimensions, net depth in in. (mm)	18 9/16 (471)
External dimensions, gross height in in. (mm)	5 15/16 (150)
External dimensions, gross width in in. (mm)	9 1/16 (230)
External dimensions, gross depth in in. (mm)	19 5/16 (490)
Net weight in lb (kg)	5 (2.30)
Gross weight in lb (kg)	7 (3.15)