

E 416  
 Insert  
 for optimum loading of 6 single and two-piece specula.



- Suitable for 6 single and two-part specula
- 7 holders 40 mm apart
- For use in upper (only with single-part specula) and lower basket

EAN: 4002516854135 / Material number: 12669750 / Old Material Number: 69541602D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8534	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Insert for gynaecology instruments	•
<b>Application area</b>	
Reprocessing of GYN instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Capacity</b>	
GYN specula [number]	6
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	6 3/16 (157)
External dimensions, net width in in. (mm)	7 1/16 (178)
External dimensions, net depth in in. (mm)	11 (279)
External dimensions, gross height in in. (mm)	8 5/16 (210)
External dimensions, gross width in in. (mm)	7 7/8 (200)
External dimensions, gross depth in in. (mm)	14 3/8 (365)
Net weight in lb (kg)	2 (0.69)
Gross weight in lb (kg)	2 (0.97)