

E 458

Insert

for the optimum loading of screw-cap teats.



- Suitable for 36 screw connection teats
- 36 compartments at 12 x 12
- equipped with a flip-up and lockable lid

EAN: 4002516854265 / Material number: 12669800 /
Old Material Number: 69545802D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other inserts	•
Application area	
Reprocessing of baby bottles	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	36
Dimensions and weight	
External dimensions, net height in in. (mm)	2 1/2 (63)
External dimensions, net width in in. (mm)	8 1/2 (215)
External dimensions, net depth in in. (mm)	17 9/16 (445)
External dimensions, gross height in in. (mm)	3 9/16 (90)
External dimensions, gross width in in. (mm)	9 1/16 (230)
External dimensions, gross depth in in. (mm)	18 5/16 (465)
Net weight in lb (kg)	3 (1.50)
Gross weight in lb (kg)	4 (1.77)