

E 466

Injector nozzle

for the optimum loading of breathing bags 8 x 333 mm.



EAN: 4002516854920 / Material number: 12669890 /
Old Material Number: 69746602D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other accessories	•
Application area	
Reprocessing anaesthetic instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net depth in in. (mm)	13 1/8 (333)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	5 15/16 (150)
External dimensions, gross depth in in. (mm)	7 7/8 (200)
Net weight in lb (kg)	2 (0.85)
Gross weight in lb (kg)	2 (1.01)