

E 494  
Insert  
for the gentle holding of 5 microtiter plates.



EAN: 4002516853947 / Material number: 12669920 /  
Old Material Number: 69549402D

<b>Product category</b>	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
SlimLine laboratory glassware washers	•
Machines for cleaning industrial parts	•
Large-chamber laboratory glasswashers for cleaning industrial parts	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 7111 [IPC]	•
PLW 8505	•
<b>Product category</b>	
Other inserts	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Capacity</b>	
Microtiter plates [number]	5
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	1 7/16 (35)
External dimensions, net width in in. (mm)	8 1/8 (205)
External dimensions, net depth in in. (mm)	17 3/8 (440)
External dimensions, gross height in in. (mm)	3 9/16 (90)
External dimensions, gross width in in. (mm)	9 1/16 (230)
External dimensions, gross depth in in. (mm)	18 5/16 (465)

E 494  
Insert  
for the gentle holding of 5 microtiter plates.



EAN: 4002516853947 / Material number: 12669920 /  
Old Material Number: 69549402D

Net weight in lb (kg)	1 (0.49)
Gross weight in lb (kg)	2 (0.72)