

A 612 Module

For the optimum reprocessing of large-volume laboratory glassware.



- Can be individually equipped with A 846, A 847, and A 848 injector nozzles
- Suitable for 2–10 l laboratory flasks, 2–6 l round flasks
- Suitable for 2–5 l Erlenmeyer flasks, 2–5 l measuring flasks
- For use in A 503 and APLW 511 mobile units

EAN: 4002516854173 / Material number: 12670000 / Old Material Number: 69661201D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert / module for narrow necked glassware	•
Module for large-volume laboratory glassware	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Capacity	
Large-volume laboratory glassware [number]	4
Dimensions and weight	
External dimensions, net height in in. (mm)	8 3/4 (221)
External dimensions, net width in in. (mm)	10 (253)
External dimensions, net depth in in. (mm)	23 3/4 (602)
External dimensions, gross height in in. (mm)	10 (253)
External dimensions, gross width in in. (mm)	11 9/16 (293)
External dimensions, gross depth in in. (mm)	25 9/16 (648)
Net weight in lb (kg)	7 (3.23)
Gross weight in lb (kg)	8 (3.75)