


A 613 Module

For the optimum reprocessing of volumetric and measuring pipettes.



- Suitable for 98 pipettes
- Suitable for measuring pipettes with a volume of 1–50 ml
- Suitable for volumetric pipettes with a volume of 1–100 ml
- For use in A 503 mobile unit
- Capacity dependent on max. diameter of pipettes

EAN: 4002516854128 / Material number: 12670020 / Old Material Number: 69661301D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert / module for pipettes	•
Application area	
Reprocessing of pipettes / phials	•
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Capacity	
Pipettes [number]	98
Dimensions and weight	
External dimensions, net height in in. (mm)	7 11/16 (194)
External dimensions, net width in in. (mm)	10 (254)
External dimensions, net depth in in. (mm)	23 13/16 (604)
External dimensions, gross height in in. (mm)	8 9/16 (217)
External dimensions, gross width in in. (mm)	11 3/16 (283)
External dimensions, gross depth in in. (mm)	24 1/4 (615)
Net weight in lb (kg)	6 (2.69)
Gross weight in lb (kg)	7 (3.06)