

E 496

Injector nozzle
with retaining clamp for intubation material.



EAN: 4002516854579 / Material number: 12670260 / Old Material Number: 69749602D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Injector nozzle	•
Application area	
Reprocessing anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Dimensions and weight	
External dimensions, net height in in. (mm)	4 3/4 (120)
External dimensions, gross height in in. (mm)	0 1/2 (12)
External dimensions, gross width in in. (mm)	3 3/16 (80)
External dimensions, gross depth in in. (mm)	4 3/4 (120)
Net weight in lb (kg)	0 (0.01)
Gross weight in lb (kg)	0 (0.02)