


E 803

Insert

for the optimum loading of up to 30 ear and nose specula.



- 160 compartments approx. 13 x 13 mm
- Mesh width on base: 1,7 mm
- For use in upper and lower basket

EAN: 4002516854715 / Material number: 12670290 / Old Material Number: 69580302D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8534	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for ENT instruments	•
Application area	
Reprocessing of ENT instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	160
Dimensions and weight	
External dimensions, net height in in. (mm)	2 1/2 (63)
External dimensions, net width in in. (mm)	6 1/2 (165)
External dimensions, net depth in in. (mm)	12 1/2 (317)
External dimensions, gross height in in. (mm)	5 1/4 (132)
External dimensions, gross width in in. (mm)	9 3/16 (233)
External dimensions, gross depth in in. (mm)	19 1/16 (484)
Net weight in lb (kg)	2 (0.80)
Gross weight in lb (kg)	2 (1.07)