

A 620

Injector module

For the optimum reprocessing of laboratory glassware, 10 injector positions.



- Can be equipped with up to 10 nozzles based on individual needs
- For use in A 503 and APLW 511 mobile units
- Part of SimpleLoad system

EAN: 4002516853909 / Material number: 12678100 / Old Material Number: 69662001D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert / module for narrow necked glassware	•
Insert for various bottles	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Part of SimpleLoad system	•
Dimensions and weight	
External dimensions, net height in in. (mm)	4 1/2 (113)
External dimensions, net width in in. (mm)	5 1/4 (133)
External dimensions, net depth in in. (mm)	23 15/16 (608)
External dimensions, gross height in in. (mm)	5 3/4 (145)
External dimensions, gross width in in. (mm)	5 15/16 (150)
External dimensions, gross depth in in. (mm)	24 7/16 (620)
Net weight in lb (kg)	4 (1.59)
Gross weight in lb (kg)	4 (1.82)