

A 621

Injector module

For the optimum reprocessing of laboratory glassware, 20 injector positions.



- Can be equipped with up to 20 nozzles based on individual needs
- For use in A 503 and APLW 511 mobile units
- Part of SimpleLoad system

EAN: 4002516855125 / Material number: 12678110 / Old Material Number: 69662101D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert / module for narrow necked glassware	•
Insert for various bottles	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Part of SimpleLoad system	•
Capacity	
Narrow necked glassware [number]	20
Dimensions and weight	
External dimensions, net height in in. (mm)	4 1/2 (113)
External dimensions, net width in in. (mm)	6 13/16 (173)
External dimensions, net depth in in. (mm)	23 15/16 (608)
External dimensions, gross height in in. (mm)	6 3/4 (170)
External dimensions, gross width in in. (mm)	8 1/4 (208)
External dimensions, gross depth in in. (mm)	25 7/16 (645)
Net weight in lb (kg)	4 (1.93)
Gross weight in lb (kg)	5 (2.27)