

E 907/1

Mesh tray

for the optimum loading of different small-scale instruments.



- With lid for small parts
- Mesh width 3 x 1 mm
- Hook for hanging in E 905/1 and E 903/1

EAN: 4002516854258 / Material number: 12678170 / Old Material Number: 69590703D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Mesh tray	•
Application area	
Reprocessing of dental instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	1 15/16 (48)
External dimensions, net width in in. (mm)	5 1/8 (129)
External dimensions, net depth in in. (mm)	6 3/4 (170)
External dimensions, gross height in in. (mm)	3 (75)
External dimensions, gross width in in. (mm)	6 5/16 (160)
External dimensions, gross depth in in. (mm)	9 5/8 (243)
Net weight in lb (kg)	1 (0.26)
Gross weight in lb (kg)	1 (0.37)