

A 622

Injector module

For the optimum reprocessing of laboratory glassware, 36 injector positions.



- Can be equipped with up to 36 nozzles
- For use in A 503 and APLW 511 mobile units
- Part of SimpleLoad system

EAN: 4002516854531 / Material number: 12678490 / Old Material Number: 69662201D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Insert / module for narrow necked glassware	•
Insert for various bottles	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Part of SimpleLoad system	•
<b>Capacity</b>	
Narrow necked glassware [number]	36
Flasks [number]	36
Round flasks [number]	36
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	4 1/2 (113)
External dimensions, net width in in. (mm)	8 1/16 (204)
External dimensions, net depth in in. (mm)	23 15/16 (608)
External dimensions, gross height in in. (mm)	5 3/4 (145)
External dimensions, gross width in in. (mm)	8 3/4 (222)
External dimensions, gross depth in in. (mm)	24 13/16 (629)
Net weight in lb (kg)	5 (2.23)
Gross weight in lb (kg)	6 (2.61)