


A 623

Module

For the optimum loading of measuring cylinders.



- Equipped with 4 nozzles
- For measuring cylinders with 1–2 l capacity
- Pivoting supports for ease of use

EAN: 4002516853893 / Material number: 12678510 / Old Material Number: 69662301D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Insert / module for measuring cylinder	•
<b>Application area</b>	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Capacity</b>	
Measuring cylinders [number]	4
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	17 1/4 (437)
External dimensions, net width in in. (mm)	10 1/16 (255)
External dimensions, net depth in in. (mm)	23 15/16 (608)
External dimensions, gross height in in. (mm)	10 5/8 (269)
External dimensions, gross width in in. (mm)	15 3/4 (400)
External dimensions, gross depth in in. (mm)	24 1/4 (615)
Net weight in lb (kg)	5 (2.42)
Gross weight in lb (kg)	8 (3.42)