

TC 1150

Upper enclosure

For 1,150 mm-wide appliances without a steam condenser or drying unit.



- Includes lockable service door
- Wall panelling up to ceiling must be provided on site

EAN: 4002516853848 / Material number: 12680670 /
Old Material Number: 69252005D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8617	•
Product category	
Cover	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	39 1/2 (1002)
External dimensions, net width in in. (mm)	32 7/8 (835)
External dimensions, net depth in in. (mm)	47 1/4 (1200)
External dimensions, gross height in in. (mm)	39 1/2 (1002)
External dimensions, gross width in in. (mm)	33 7/8 (860)
External dimensions, gross depth in in. (mm)	47 1/4 (1200)
Net weight in lb (kg)	133 (60.13)
Gross weight in lb (kg)	166 (75.13)