

MC/2

Transfer cart

height adjustable, for loading & unloading washer-&disinfectors.



- Equipped with 4 castors, 2 of which can be locked
- 2 loading racks (trough shaped) as well as handle and cover plate

EAN: 4002516853855 / Material number: 12680970 / Old Material Number: 69200113D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Storage / transport system	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	39 1/2 (1003)
External dimensions, net width in in. (mm)	23 5/8 (600)
External dimensions, net depth in in. (mm)	31 3/4 (805)
External dimensions, gross height in in. (mm)	47 11/16 (1210)
External dimensions, gross width in in. (mm)	31 1/2 (800)
External dimensions, gross depth in in. (mm)	47 1/4 (1200)
Net weight in lb (kg)	40 (18.00)
Gross weight in lb (kg)	105 (47.80)