

MT
Trolley
for loading & unloading washer-&disinfectors.

EAN: 4002516854159 / Material number: 12681000 / Old Material Number: 69200114D



- For the storage and transport of baskets and inserts
- Divided into 4 levels that can be adjusted by 102.5 mm
- Moveable with 4 lockable castors
- Loading volume 549 x 599 mm

Product category	
Glassware Washers, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Washer-disinfectors, fire brigades, and emergency services	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWA 8672	•
PWD 8682	•
PWD 8682 C	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Storage / transport system	•
Dimensions and weight	
External dimensions, net height in in. (mm)	78 3/16 (1985)
External dimensions, net width in in. (mm)	26 1/8 (662)
External dimensions, net depth in in. (mm)	24 5/16 (616)
External dimensions, gross height in in. (mm)	36 1/4 (920)
External dimensions, gross width in in. (mm)	30 9/16 (775)
External dimensions, gross depth in in. (mm)	78 13/16 (2001)
Net weight in lb (kg)	51 (23.00)
Gross weight in lb (kg)	83 (37.60)