

A 306/1 Module

For the optimum loading of large-volume laboratory glassware.



- Suitable for 4 1–2 litre measuring cylinders, high shape
- Suitable for 2 1–2 litre measuring cyl., low shape + 2 high shape
- Plastic-coated contact surfaces
- Measuring cylinders even removable when module is in machine
- Maximum length of the measuring cylinders in the A 150 lower basket: 500 mm
- Up to 2 modules can be inserted in baskets A 150 and APLW 106

EAN: 4002516854456 / Material number: 12683100 /
Old Material Number: 69630602D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Insert / module for measuring cylinder	•
Module for large-volume laboratory glassware	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Product characteristics	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Colour	Stainless Steel, Black, White
Capacity	
Measuring cylinders [number]	4
Dimensions and weight	
External dimensions, net height in in. (mm)	9 5/16 (235)
External dimensions, net width in in. (mm)	16 1/2 (418)
External dimensions, net depth in in. (mm)	18 9/16 (471)
External dimensions, gross height in in. (mm)	9 7/8 (250)
External dimensions, gross width in in. (mm)	15 3/8 (390)
External dimensions, gross depth in in. (mm)	19 1/8 (485)
Net weight in lb (kg)	6 (2.88)
Gross weight in lb (kg)	8 (3.42)