

A 501

Mobile unit

With 4 levels for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays
- Extendable to 5 or 6 levels with A 601 racks

EAN: 4002516855361 / Material number: 12684440 / Old Material Number: 69650102D

| | |
|--|--------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| Product category | |
| Mobile unit for inserts | • |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black |
| Capacity | |
| DIN mesh trays [number] | 18 |
| SPRI mesh trays [number] | 12 |
| Dimensions and weight | |
| External dimensions, net height in in. (mm) | 26 (659) |
| External dimensions, net width in in. (mm) | 25 5/8 (650) |
| External dimensions, net depth in in. (mm) | 28 (710) |
| External dimensions, gross height in in. (mm) | 28 (710) |
| External dimensions, gross width in in. (mm) | 31 1/2 (800) |
| External dimensions, gross depth in in. (mm) | 33 7/8 (860) |
| Net weight in lb (kg) | 46 (21.00) |
| Gross weight in lb (kg) | 72 (32.80) |