

A 503

Mobile unit

No levels, for cleaning laboratory glassware/utensils.



- For loading modules and inserts
- Depending on the height, up to 6 modules can be used
- Upper rail and upper strut can be removed

EAN: 4002516854784 / Material number: 12685380 /
Old Material Number: 69650302D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Mobile unit with module connection	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Dimensions and weight	
External dimensions, net height in in. (mm)	17 1/16 (433)
External dimensions, net width in in. (mm)	31 1/8 (790)
External dimensions, net depth in in. (mm)	25 5/8 (650)
External dimensions, gross height in in. (mm)	31 1/8 (790)
External dimensions, gross width in in. (mm)	28 3/16 (715)
External dimensions, gross depth in in. (mm)	33 7/8 (860)
Net weight in lb (kg)	25 (11.40)
Gross weight in lb (kg)	49 (22.20)