

E 491
Insert
for the optimum loading of rotating instruments.



- Suitable for 18 turbine shafts
- Suitable for 30 contra angle shafts
- Also suitable for rotating instruments, such as used in podiatry

EAN: 4002516855194 / Material number: 12686260 / Old Material Number: 68549103D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Large-capacity washer-disinfectors, dental	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
PWA 8672	•
PWD 7121	•
PWD 7122	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert / module for dental instruments	•
Application area	
Reprocessing of dental instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Capacity	
Sheaths [number]	48
Dimensions and weight	
External dimensions, net height in in. (mm)	2 1/8 (53)
External dimensions, net width in in. (mm)	2 13/16 (70)
External dimensions, net depth in in. (mm)	4 3/4 (120)
External dimensions, gross height in in. (mm)	3 15/16 (100)
External dimensions, gross width in in. (mm)	10 7/16 (265)
External dimensions, gross depth in in. (mm)	7 1/8 (180)
Net weight in lb (kg)	1 (0.30)
Gross weight in lb (kg)	2 (1.00)