

APLW 909
 Installation kit
 for dispensing neutralizing agents from 10-litre media canisters.



- For using 10-litre media canisters outside the machine
- Suction lance for neutralizing agents

EAN: 4002516930044 / Material number: 12755470 / Old Material Number: 69242021D

Product category	
SlimLine laboratory glassware washers	•
SlimLine laboratory glasswashers for cleaning industrial parts	•
Machine affiliation	
PLW 7111	•
PLW 7111 [Safety]	•
PWD 7121	•
PWD 7122	•
Product category	
Conversion kit	•
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of dental instruments	•
Product characteristics	
Material	Plastic
Dimensions and weight	
External dimensions, gross height in in. (mm)	3 3/4 (95)
External dimensions, gross width in in. (mm)	6 3/4 (170)
External dimensions, gross depth in in. (mm)	21 11/16 (550)
Net weight in lb (kg)	1 (0.34)
Gross weight in lb (kg)	1 (0.50)