

ProCare Dent 10 A - 5 l
Liquid detergent, alkaline, 5 l
for the automatic reprocessing of dental instruments.



- Suitable for dental instruments, not for handpieces
- Adapted to the process sequences in Miele washer-disinfectors
- Excellent material/corrosion protection for stainless steel instruments
- Does not contain surfactants
- Free from oxidative additives

EAN: 4029577158194 / Material number: 12786490 / Old Material Number: 66912201V6

Product category	
Thermal disinfectors	•
Applications	
Removal of organic residues	•
Dispensing	
pH value at 20°C (68°F)	14,0
Dispensing recommendation	
Dispensing recommendation per rinse	3 ml/l (0.3%)
Dispensing recommendation (per cycle)	3 ml/l (0.3%)
Large opening for suction lance	•
Intensity of soiling	Strong
Suitable for all levels of water hardness	•
Main applications	
Dental practices	•
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Phosphate, Silicates
Minimum storage temperature in °C	-15
Maximum storage temperature in °C	30
Dimensions and weight	
External dimensions, gross depth in mm	150.00
Product characteristics	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Phosphate, Silicates
Minimum storage temperature in °C	-15
Maximum storage temperature in °C	30
Suitable for standard instruments	•
Suitable for automatic reprocessing of instruments and utensils	•
Does not contain oxidants	•
High degree of material protection against corrosion	•
Good cleaning capability	•
Does not contain surfactants	•
Not suitable for handpieces	Suitable for standard instruments, Not suitable for transmission instr., For machine reprocessing of instruments
Dimensions and weight	
External dimensions, net height in in. (mm)	8 9/16 (217)
External dimensions, net width in in. (mm)	5 5/8 (142)
External dimensions, net depth in in. (mm)	9 13/16 (248)
External dimensions, gross height in in. (mm)	9 7/8 (250)
External dimensions, gross width in in. (mm)	7 1/2 (190)
External dimensions, gross depth in in. (mm)	5 15/16 (150)
Gross weight in lb (kg)	0.35 (10)

ProCare Dent 10 A - 5 l
Liquid detergent, alkaline, 5 l
for the automatic reprocessing of dental instruments.



EAN: 4029577158194 / Material number: 12786490 / Old Material Number: 66912201V6

Fill volume in fl.oz. (ml)	169.07 (5000) fl oz (US) (Milliliter)
----------------------------	--

Compliance, test certificates, and symbols	
CE	•
GHS / CLP compliant	•